

Treatment	Replicate	Initial Log₁₀ Inoculated	Final Log₁₀ Recovered	Log₁₀ Reduction
Ozone	1	8.20	1.78	6.43
	2	8.20	1.78	6.43
	3	8.20	1.48	6.73

Table 2. Inactivation of P22 bacteriophages on N95 mask coupons using ozone.

Treatment	Replicate	Initial Log₁₀ Inoculated	Final Log₁₀ Recovered	Log₁₀ Reduction
Vaporized Hydrogen Peroxide	1	8.20	1.48	6.73
	2	8.20	1.78	6.43
	3	8.20	1.48	6.73

Table 2. Inactivation of P22 bacteriophages on N95 mask coupons using VHP sterilization system.